



A.T.P RING
CAPSULAR TENSION RING

Description
One piece PMMA ring
Ring featuring new in special design
ETO sterilized



ATP CTRs are used to stabilize and recenter the capsular bag in the following cases:

- Significant anterior weakness
- Lens subluxation
- Pseudophakia
- Myopia
- Cataracts
- Soft IOL
- Marfan's Syndrome

Specification

Material: Blue PMMA/Clear PMMA

Model	Size	Thickness
CTR-12	10-12	0.20 mm
CTR-13	11-13	0.30 mm

