



A.T.P. INJECTOR

ATP Injection system™ Type A & C

These lens delivery systems are designed for 1-piece lenses with incision size of <3.0 mm and <2.0 mm corneals:

- One single-use injector
- Cartridge FlyGlide™
- Silicon stopper

There are two Injector Systems for the different IOL types

A.T.P. injector Systems make eye surgery safer, simpler and more effective. The design is highly efficient, cost-effective delivery system and it is achieved by Consider customer's needs and expectation, especially about current new MICS technique.

The surgeon can insert Foldable Lens with different incision size with various A.T.P. injectors.

